

Archwires, Straight Lengths and Springs



A variation of Nitinol Archwires
A variation of Stainless Steel Archwires
Several Straight Lengths
A large variety of Nitinol Springs

Heat-Activated

This wire is heat activated by mouth temperature. It delivers light continuous forces, moving teeth with less discomfort to patient.

Arch form Medium



Round	Upper	Lower
.014	536-214	536-314
.016	536-216	536-316
.018	536-218	536-318
.020	536-220	536-320



Rectangular	Upper	Lower
.016 x .016	536-261	536-361
.016 x .022	536-262	536-362
.017 x .025	536-264	536-364
.018 x .025	536-266	536-366
.019 x .025	536-267	536-367
.021 x .025	536-268	536-368

Reverse Curve

This is a malleable and easy to install wire. It is used mainly for the correction of accentuated curve of spee.

Arch form Medium



Round	Upper	Lower
.012	517-812	517-912
.014	517-814	517-914
.016	517-816	517-916
.018	517-818	517-918
.020	517-820	517-920



Rectangular	Upper	Lower
.016 x .016	517-860	517-960
.016 x .022	517-862	517-962
.017 x .025	517-864	517-964
.018 x .025	517-865	517-965
.019 x .025	517-867	517-967
.021 x .025	517-868	517-968

Super-Elastic

It is easy to engage and provides light continuous forces. It can be used after engaging Heat-Activated wires. This wire allows good torque control.

Arch form Medium



Round	Upper	Lower
.012	510-212	510-312
.014	510-214	510-314
.016	510-216	510-316
.018	510-218	510-318
.020	510-220	510-320



Rectangular	Upper	Lower
.016 x .016	510-261	510-361
.016 x .022	510-262	510-362
.017 x .025	510-264	510-364
.018 x .025	510-266	510-366
.019 x .025	510-267	510-367
.021 x .025	510-268	510-368

Arch form Large



Round	Upper	Lower
.012	517-212	517-312
.014	517-214	517-314
.016	517-216	517-316
.018	517-218	517-318
.020	517-220	517-320



Rectangular	Upper	Lower
.016 x .016	517-261	517-361
.016 x .022	517-262	517-362
.017 x .025	517-264	517-364
.018 x .025	517-266	517-366
.019 x .025	517-267	517-367
.021 x .025	517-268	517-368

Stainless Steel

Manufactured with high quality 304-V steel. This special metal provides optimum flexibility, decreasing breakage.

Arch form Medium



Round	Upper	Lower
.012	522-812	522-912
.014	522-814	522-914
.016	522-816	522-916
.018	522-818	522-918
.020	522-820	522-920



Rectangular	Upper	Lower
.016 x .016	522-861	522-961
.016 x .022	522-862	522-962
.017 x .025	522-864	522-964
.018 x .025	522-866	522-966
.019 x .025	522-867	522-967
.021 x .025	522-868	522-968

Arch form Large



Round	Upper	Lower
.012	522-212	522-312
.014	522-214	522-314
.016	522-216	522-316
.018	522-218	522-318
.020	522-220	522-320



Rectangular	Upper	Lower
.016 x .016	522-261	522-361
.016 x .022	522-262	522-362
.017 x .025	522-264	522-364
.018 x .025	522-266	522-366
.019 x .025	522-267	522-367
.021 x .025	522-268	522-368

Straight Lengths

Manufactured with high quality 304-V steel. This special metal provides optimum flexibility, decreasing breakage.

20 per tube



14" - Round

- .012
- .014
- .016
- .018
- .020

Code

- 502-112
- 502-114
- 502-116
- 502-118
- 502-120

10 per tube



14" - Rectangular


- .016 x .016
- .016 x .022
- .017 x .025
- .018 x .025
- .019 x .025
- .021 x .025

Code

- 502-161
- 502-162
- 502-164
- 502-166
- 502-167
- 502-168

Elgiloy - Blue Temper
Appropriate for tempering when working on the wire.

10 per tube



Rectangular

- .016 x .016
- .016 x .022

Code

- 520-161
- 520-162

Nitinol

Open Coil Spring



Force	Diameter	Length	Package	Code
Average	010x030	17 cm	3 per tube	410-142

Distalizing Spring - Open and Closed



Force	Diameter	Length	Package	Code
Average	010x030	17 cm	3 per tube	410-100

Coil Spring



Force	Diameter	Length	Package	Code
Variable	010x030	9 mm	10 per box	410-108
Variable	010x030	12 mm	10 per box	410-109

Elastics



'O' Ligatures
Separator Modules
Elastomeric Ligatures
Elastomeric Chains
Intraoral Elastics



'O' Ligature

Diameter .110








1.000 ligatures per pack

*Colors depicted are similar to the actual product.

Color*

Code


	Yellow	450-110
	Baby Blue	450-128
	Light Blue	450-115
	Dark Blue	450-105
	Navy Blue	450-126
	Gray	450-102
	Clear	450-101
	Assorted Package	450-100
	Orange	450-109
	Lilac	450-112
	Pink	450-106
	Silver	450-116
	Black	450-108
	Rose	450-127
	Purple	450-107
	Dark Green	450-111
	Metallic Green	450-119
	Extreme Green	450-120
	Red	450-104
	Violet	450-118

Separator Modules

Diameter

Color*

Code





960 separators per pack

*Colors depicted are similar to the actual product.

	Blue	466-424
	Blue	466-426 (Radiopaque)

Elastomeric Ligatures

Elastomeric ligatures are available with two options: patient stick and universal stick. Both have a .115 diameter and each package contains, approximately, 1000 elastomeric ligatures.

Color*	Patient Stick	Universal Stick
	 <ul style="list-style-type: none"> • 24 elastomeric ligatures • 42 sticks per pack • 1008 elastomeric ligatures per pack 	 <ul style="list-style-type: none"> • 12 sticks per pack • 1008 elastomeric chains per pack
Yellow	466-610	466-710
Baby Blue	466-628	466-728
Light Blue	466-615	466-715
Dark Blue	466-605	466-705
Dark Blue Metallic	466-614	466-714
Mediterranean Blue	466-629	466-729
Navy Blue	466-626	466-726
White	466-603	466-703
Gray	466-602	466-702
Clear	466-601	466-701
Assorted Package	466-600	466-700
Orange	466-609	466-709
Lilac	466-612	466-712
Pink	466-606	466-706
Silver	466-616	466-716
Black	466-608	466-708
Rose	466-627	466-727
Baby Rose	466-624	466-724
Fuschia	466-625	466-725
Purple	466-607	466-707
Green	466-621	466-721
Light Green	466-622	466-722
Metallic Green	466-619	466-719
Dark Green	466-611	466-711
Extreme Green	466-620	466-720
Red	466-604	466-704
Red Metallic	466-623	466-723
Violet	466-618	466-718

*Colors depicted are similar to the actual product.

Elastomeric Chains

Elastomeric chains are available in 10 different colors and with three different sizes: small, medium and large. They can be placed in our alastic dispenser.

4.6 meters per bobbin



Color*	 Small	 Medium	 Large
 Yellow	467-710	467-810	467-910
 Light Blue	467-715	467-815	467-915
 Dark Blue Metallic	467-714	467-814	467-914
 Dark Blue	467-705	467-805	467-905
 Navy Blue	467-726	467-826	467-926
 White	467-703	467-803	467-903
 Gray	467-702	467-802	467-902
 Clear	467-701	467-801	467-901
 Golden	467-717	467-817	467-917
 Orange	467-709	467-809	467-909
 Lilac	467-712	467-812	467-912
 Pink	467-706	467-806	467-906
 Silver	467-716	467-816	467-916
 Black	467-708	467-808	467-908
 Rose	467-725	467-825	467-925
 Purple	467-707	467-807	467-907
 Green	467-721	467-821	467-921
 Light Green	467-722	467-822	467-922
 Dark Green	467-711	467-811	467-911
 Metallic Green	467-719	467-819	467-919
 Emerald Green	467-723	467-823	467-923
 Extreme Green	467-720	467-820	467-920
 Red	467-704	467-804	467-904
 Violet	467-718	467-818	467-918

*Colors depicted are similar to the actual product.

Intraoral

They provide stable and continuous force.



Diameter	LIGHT Force 2.5oz = 71g	MEDIUM Force 4.5oz = 127g	HEAVY Force 6.5oz = 184g	Package
1/8	464-100	464-102	464-104	1000 per pack
3/16	464-200	464-202	464-204	1000 per pack
1/4	464-300	464-302	464-304	1000 per pack
5/16	464-400	464-402	464-404	1000 per pack
3/8	464-500	464-502	464-504	1000 per pack
1/2	463-200	463-202	463-204	100 per pack



Attachments




Palatal Bar Sheaths
Spur
Lingual Button
Ram Cleats Hooks



Palatal Bar Sheaths

	Width	Code
 <p>Palatal Bar Lingual Sheath with window and hook.</p>	.036 - U/R, L/L	425-002
 <p>Palatal Bar Lingual Sheath with window and lock, without hook.</p>	.036 - U/L, L/R .036 - U/R, U/L, L/L, LR .036 - U/R, L/L	425-003 425-011 425-013

 <p>Palatal Bar Lingual Sheath with window and lock, with hook (mesial or distal).</p>	.036 - U/L, L/R	425-014
--	-----------------	---------

Spur




	Package	Code
	10 per pack	425-810



Nogueira® Bondable Spur. A new proposal when treating atypical deglutition.
Dr. Francisco Ferreira Nogueira


Lingual Button

Lingual button mushroom head highly polished with suitable profile to hold elastics.

Bondable Lingual Buttons 10 per pack	Base	Code
	.036 - U/L, L/R	400-201
	.036 - U/R, U/L, L/L, LR	400-202
Welding Lingual Buttons 10 per pack		
	Flat	400-005
	Curved	400-006
	Over curved	400-007
Bondable mini-buttons 10 per pack	Package	Code
	10 per pack	400-213

Ram Cleats Hooks

Versatile lingual attachments that can be used as a hook or as a seating cleat.

Welding 10 per pack	Position	Code
	Anterior, Cuspid and Bicuspid	400-002
	Molar	400-003